<b>DA</b>	<b>Pædiatrisk</b>
<b>Genanvendeligt hyper-hypotermisk tæppe</b>	<b>Herbruikbare hyper-/hypothermiedeken</b>
<b>Termisk reguleringsystem</b>	<b>Thermische regelsysteem</b>

**ADVARSLER**

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<b>CSZ</b>	<b>REF</b>	<b>194CPC</b>
<b>A GENTHERM COMPANY</b>	<b>Size:</b> 22in x 30in (55.9cm x 76.2cm)	

# PlastiPad® Pediatric Reusable Hyper-Hypothermia Blanket Thermal Regulating System

**WARNINGS**

- Licensed healthcare practitioner’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominences. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and PlastiPad®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the PlastiPad® blanket as punctures may occur. Do not use pins to secure Plastipad® blanket.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.

**PRECAUTIONS**

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the blanket while the blanket is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The PlastiPad® has been designed and tested for use with CSZ thermal regulating systems only.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients’ skin during application of prolonged heat.
- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

<b>CSZ</b>	<b>EC REP</b>	Authorized European Representative
<b>A GENTHERM COMPANY</b>	<b>C</b>	CEpartner4U

<b>CSZ</b>	<b>0344</b>	Authorized European Representative
<b>A GENTHERM COMPANY</b>	<b>C</b>	CEpartner4U

Cincinnati Sub-Zero Products, LLC  
12011 Mosteller Road  
Cincinnati, OH 45241  
Tel: 1-800-989-7373  
[www.cszmedical.com](http://www.cszmedical.com)

<b>CSZ</b>	<b>EC REP</b>	Authorized European Representative
<b>A GENTHERM COMPANY</b>	<b>C</b>	CEpartner4U

Esdorndlaan 13  
3951 DB Maarn  
The Netherlands  
[www.CEpartner4U.com](http://www.CEpartner4U.com)

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