DA	Pædiatrisk
Genanvendeligt hyper-hypotermisk tæppe	Herbruikbare hyper-/hypothermiedeken
Termisk reguleringssystem	Thermische regelsysteem

ADVARSLER

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
CSZ	CSZ
A GENTHERM COMPANY	A GENTHERM COMPANY

REF	196CPC
Size: 24in x 60in (61.0cm x 152.4cm)	

PlastiPad®	PlastiPad®
Adult	Adult
Reusable Hyper-Hypothermia Blanket	Reusable Hyper-Hypothermia Blanket
Thermal Regulating System	Thermal Regulating System

WARNINGS

- Licensed healthcare practitioner’s order is required for temperature setting and for use of the therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominences. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and PlastiPad®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the PlastiPad® blanket as punctures may occur. Do not use pins to secure Plastipad® blanket.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.

PRECAUTIONS

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the blanket while the blanket is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The PlastiPad® has been designed and tested for use with CSZ thermal regulating systems only.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients’ skin during application of prolonged heat.
- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

CSZ	CSZ
A GENTHERM COMPANY	A GENTHERM COMPANY

EC REP	EC REP
Authorized European Representative	Authorized European Representative

CE	CE
0344	0344

CSZ	CSZ
A GENTHERM COMPANY	A GENTHERM COMPANY
Cincinnati Sub-Zero Products, LLC	Cincinnati Sub-Zero Products, LLC
12011 Mosteller Road	12011 Mosteller Road
Cincinnati, OH 45241	Cincinnati, OH 45241
Tel: 1-800-989-7373	Tel: 1-800-989-7373
www.cszmedical.com	www.cszmedical.com

CSZ	CSZ
A GENTHERM COMPANY	A GENTHERM COMPANY
Cincinnati Sub-Zero Products, LLC	Cincinnati Sub-Zero Products, LLC
12011 Mosteller Road	12011 Mosteller Road
Cincinnati, OH 45241	Cincinnati, OH 45241
Tel: 1-800-989-7373	Tel: 1-800-989-7373
www.cszmedical.com	www.cszmedical.com

CE	CE
0344	0344

Authorized European Representative	Authorized European Representative
CEpartner4U	CEpartner4U
Esdoornlaan 13	Esdoornlaan 13
3951 DB Maarn	3951 DB Maarn
The Netherlands	The Netherlands
www.CEpartner4U.com	www.CEpartner4U.com

57562-C	57562-C
(09/2018)	(09/2018)

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