Voksen

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Thermal Regulating System

WARNINGS

- Physicians order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a physician. Higher risk patients; including pediatric, temperature sensitive patients and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to all bony prominences. Notify physician for changes in patient condition.
- Follow physician order/hospital policy for discontinuation of therapy. Notify physician for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and PlastiPad®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the PlastiPad® blanket as punctures may occur. Do not use pins to secure PlastiPad® blanket.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.

PRECAUTIONS

- Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare professional.
- Presence of water or gel on the blanket while the blanket is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The PlastiPad® has been designed and tested for use with CSZ thermal regulating systems only.

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Size: 24in x 60in
(61.0cm x 152.4cm)

INDICATIONS FOR USE

The PlastiPad® reusable hyper-hypothermia blanket is intended to be used by trained healthcare professionals as part of CSZ's thermal regulating systems, specifically the Blanketrol® II, Blanketrol® III, Hemotherm®, or Norm-O-Temp® to lower or to raise a patient's temperature and/or maintain a desired patient temperature through conductive heat transfer.

CONTRAINDICATIONS

None known.

SET UP

- Unfold blanket on bed, stretcher, or operating room (OR) table with hose routed without kinks towards the unit. Check blanket before each use for product integrity.
- Check the equipment reservoir to ensure that it is full of distilled water.
- Connect male coupling and female coupling from the Cat#286 hose to unit. Make sure they "click" as they lock together.
- Attach blanket to the Cat#286 hose by inserting the blanket connectors into the hose connectors. Connectors should "click" as they lock together.
- Turn unit ON and select desired temperature setting.
- Place a single sheet between the patient's skin and blanket. The blanket may be filled either before or after the patient is placed on blanket. Blanket may be under or over the patient as desired.
- After filling and positioning, make sure blanket and hose are free of kinks that may restrict water flow.

OPERATING INSTRUCTIONS

- Follow the CSZ thermal regulating system operation manual for operating instructions.
- To discontinue therapy: turn the regulating system unit off. The water in the blanket will automatically drain back into the unit. Allow it to fully drain before disconnecting the blanket from unit or removing from bed, stretcher, or table.

CLEANING

- Clean blanket with a mild soap and water solution. Disinfect according to hospital/institution protocol. CSZ does not recommend sterilization. Do not autoclave.

STORAGE

- To store blanket, gently fold the corner containing the tubing towards the center of the blanket. Roll or fold blanket loosely avoiding sharp creases.
- Store blanket in cool, dry environment, free from items that may puncture blanket.

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